

### ***INSTRUMENTATION CHANGE***

On July 26th, 2011 we began using the Siemens Vista Analyzer for our chemistry testing. This change required updating the reference range for some of the tests done on the new instrument. The table below lists the tests involved and the new reference range versus the old reference range.

<b>Test</b>	<b>New Reference Range</b>		<b>Current Reference Range</b>	
	<b>Male</b>	<b>Female</b>	<b>Male</b>	<b>Female</b>
<b>ALT</b>	10-65	10-65	20 - 65	20 - 65
<b>Amylase</b>	25-115	25-115	28 - 100	28 - 100
<b>Bilirubin</b>	0.2 - 1.3	0.2 - 1.0	0.2 - 1.5	0.2 - 1.2
<b>Calcium</b>	8.5 - 10.1	8.5 - 10.1	8.5 - 10.4	8.5 - 10.4
<b>Lipase</b>	70 - 300	70 - 300	65 - 230	65 - 230
<b>Total Protein</b>	6.2 - 8.3	6.5 - 8.3	6.0 - 8.3	6.0 - 8.3
<b>Troponin I</b>	<0.05	<0.05	<0.07	<0.07